

BestPracticesSTBiocBook

Table of contents

Welcome	3
Docker image	4
RStudio Server	5
Session info	6

Welcome

Package: BestPracticesSTBiocBook **Authors:** First Last [aut, cre] **Compiled:** 2024-09-25
Package version: 0.98.0 **R version:** R version 4.4.1 (2024-06-14) **BioC version:** 3.20
License: MIT + file LICENSE

This is the website for the online book ‘**Best Practices for Spatial Transcriptomics Analysis with Bioconductor**’.

This book provides discussion and interactive examples on best practices for computational analysis workflows for spatial transcriptomics data, using the [Bioconductor](#) framework within R. The chapters contain details on individual analysis steps as well as complete example workflows, with interactive example datasets and R code.

The book is organized into several parts, including introductory materials, analysis steps, and example workflows.

Additional details on analysis workflows for non-spatial single-cell data as well as further introductory materials on R and Bioconductor can be found in the related book [Orchestrating Single-Cell Analysis with Bioconductor \(OSCA\)](#).

Docker image

A Docker image built from this repository is available here:

ghcr.io/lmweber/bestpracticesstbiocbook

Get started now

You can get access to all the packages used in this book in < 1 minute, using this command in a terminal:

Listing 0.1 bash

```
docker run -it ghcr.io/lmweber/bestpracticesstbiocbook:devel R
```

RStudio Server

An RStudio Server instance can be initiated from the **Docker** image as follows:

Listing 0.2 bash

```
docker run \  
  --volume <local_folder>:<destination_folder> \  
  -e PASSWORD=OHCA \  
  -p 8787:8787 \  
  ghcr.io/lmweber/bestpracticesstbiocbook:devel
```

The initiated RStudio Server instance will be available at <https://localhost:8787>.

Session info

 Click to expand

